

Product Brochure

Smik



For a brighter day

here's a little bit about us

our purpose

To deliver great service to our customers and improve the quality of life and promote independence to those who use our products – for a brighter day!

our history

Smik is our retail packaged home healthcare product range of Care Quip. Our team have been helping Australians live healthier, more comfortable and more independent lives since we started selling walking sticks out of a garage in Western Sydney in 1981.

Since then we have become Australia's most experienced and trusted supplier of mobility and homecare equipment to pharmacies and homecare retailers. Re-imagining the homecare offering, Smik has become a vibrant brand and offers a strong yet diverse product range of retail packaged products offered exclusively to our customers.

our people

We know we are only as good as our people – and that's why we have the best people in the business. Care is in the DNA of everything we do, and it guides every member of our team across Australia. We all share the same goal of promoting independence and improving the quality of lives through consistent, quality service and strong long-term relationships. We are dedicated to ensuring a brighter day for you and your clients!

Contents

Mobility

03



Wheelchairs

08



Bathroom

09



Daily Living

12



Pressure Care

17



Mobility

Walking Stick

- Lightweight aluminium
- Adjustable height
- Polished timber handle
- Tip size 19mm

Code	Colour	Height Approx.	Max User Weight
HH0300	Black	770 - 1000mm	100kg
HH0310	Bronze	770 - 1000mm	100kg



Patterned Walking Stick

- Lightweight aluminium
- Adjustable height
- Polished timber handle
- Tip size 19mm

Code	Style	Height Approx.	Max User Weight
HH0510	Green Patterned	770 - 1000mm	100kg
HH0520	Pink Floral	770 - 1000mm	100kg
HH0550	Blue Floral	770 - 1000mm	100kg
HH0620	Purple Floral	770 - 1000mm	100kg



Folding Walking Stick

- Lightweight aluminium
- Folds for easy storage and transportation
- Adjustable height
- Polished timber handle
- Tip size 19mm



Code	Colour	Height Approx.	Max User Weight
HH0140	Black	820 - 920mm	100kg
HH0560	Green Patterned	820 - 920mm	100kg
HH0570	Pink Floral	820 - 920mm	100kg
HH0580	Purple Floral	820 - 920mm	100kg
HH0610	Blue Floral	820 - 920mm	100kg
HH0150	Bronze	820 - 920mm	100kg

Walking Stick Strap

- Nylon woven wrist strap attaches to neck of walking stick
- Prevents the stick from being accidentally dropped
- Allows the user to freely utilise both hands while keeping the stick close and accessible
- Elastic loop fits most walking sticks

Code	Walking Stick Thickness	Colour
HZ0360	16 - 19mm	Brown
HZ0370	22 - 25mm	Black



Frog Walking Stick Holder - HZ3060

- Spring loaded plastic holder clamps onto walking stick
- Enables the walking stick to be held upright on table edge
- Foam pad assists surface grip
- Easily fitted to and removed from most walking sticks



Walking Stick Tip - HZ0006

- Patterned base for additional grip
- Supplied individually in retail packaging
- Size 19mm
- Black colour



Self Standing Cane Tip

- Suits 19mm and 22mm stick
- 6 point design
- Allows walking stick to stand upright by itself
- Adds stability for users



Code	Size	Colour
HZ3070	19mm	Black
HZ3080	22mm	Black

Quad Stick

- Adjustable height
- Foam handgrip
- Handle and stem rotate for left or right hand use
- Low profile offset base
- Lightweight aluminium stem and steel base
- Tip size 13mm



Code	Colour	Height	Base Width	Base Length	Max User Weight
HE0090	Bronze	715 - 940mm	150mm	200mm	100kg
HE0100	Silver	715 - 940mm	150mm	200mm	100kg

Underarm Crutches

- Padded top, padded handgrip and rubber tip
- Adjustable handle and leg height
- Lightweight aluminium construction
- Tip size 22mm
- Supplied as a pair in retail packaging



Code	Size	Crutch Height	Max User Weight	Approx. Patient Height
HA0210	Small	945 - 1150mm	130kg	1370 - 1570mm
HA0200	Medium	1135 - 1340mm	130kg	1570 - 1780mm
HA0190	Large	1345 - 1550mm	130kg	1780 - 1980mm

Elbow Crutches

- Adjustable cuff and leg height
- Moulded plastic handgrip
- Lightweight aluminium
- Tip size 22mm
- Supplied as a pair



Code	Size	Handgrip Height	Cuff to Handgrip Height	Max User Weight
HA0121	Medium	555 - 810mm	180 - 255mm	180kg
HA0111	Large	685 - 935mm	200 - 275mm	180kg
HA0101	Extra Large	840 - 1095mm	200 - 275mm	180kg

Crutch Tips - HZ5020

- Deep profile base for heavy use
- Patterned base for additional grip
- Supplied as a pair
- Tip size 22mm
- Grey colour



Padded Crutch Top - HZ5050

- Moulded rubber underarm crutch tops
- Supplied as pair
- To suit aluminium underarm crutches
- Grey colour



Padded Crutch Handgrip - HZ5060

- Moulded rubber crutch handle grips
- Supplied as pair
- To suit aluminium underarm crutches
- Grey colour



Triumph Walker

- Compact folding design for easy storage and transportation
- Height adjustable frame and handles
- Removable backrest
- Padded seat
- Lock on brake feature
- Vinyl storage bag



Code	Color	Seat Width	Seat Depth	Seat Height	Overall Width	Push Handle Height	Max User Weight
HF0730	Burgundy	350mm	310mm	520 - 560mm	560mm	740 - 850mm	120kg
HF0740	Blue	350mm	310mm	520 - 560mm	560mm	740 - 850mm	120kg
HF0750	Silver	350mm	310mm	520 - 560mm	560mm	740 - 850mm	120kg

Wheelchairs

Triton Self Propelled

- Swingaway removable legrests and removable armrests
- Breathable padded nylon upholstery
- Rear backrest storage pocket
- Heel loops
- Solid tyres
- Folding steel frame



Code	Color	Seat Width	Seat Depth	Seat Height	Back Height	Arm Height	Push Handle Height	Overall Width	Gross Weight	Net Weight	Max User Weight
NC0280	Burgundy	460mm	400mm	490mm	430mm	245mm	900mm	650mm	19.1kg	15kg	120kg
NC0290	Titanium	460mm	400mm	490mm	430mm	245mm	900mm	650mm	19.1kg	15kg	120kg

Vito Transit - NC1660

- Swingaway removable legrests and fixed armrests
- Desk armrests allow wheelchair to be pushed close to table
- Adjustable nylon buckle seat belt for additional patient safety
- Breathable nylon upholstery
- Rear backrest storage pocket
- Folding backrest and frame for easy storage and transportation
- Lever brake and attendant hand brakes provide more control for the carer
- Heel loops
- Solid tyres
- Lightweight aluminium frame
- Black paint finish



Seat Width	Seat Depth	Seat Height	Back Height	Arm Height	Push Handle Height	Overall Width	Gross Weight	Net Weight	Max User Weight
460mm	400mm	480mm	440mm	210mm	910mm	630mm	11.6kg	10.3kg	100kg

Bathroom

Bathboard - AA0110

- Wider contour at one end to assist transfers
- Soap dish and removable grab rail
- Integrated slots drain water away from seating area
- Adjustable feet to fit most baths



Width	Depth	Max User Weight
680mm	350mm	150kg

Shower Chair - AA0250

- Aluminium rust resistant frame
- Removable backrest
- Adjustable height
- Seat with drain holes
- Padded armrests



Seat Width	Seat Depth	Seat Height - Low Position	Seat Height - High Position	Backrest Height	Max User Weight
460mm	330mm	400mm	500mm	320mm	114kg

Bath & Shower Mats

- Provides a non-slip area within the bath or shower
- Moulded suction cups on underside for added grip
- Beige colour
- Made from quality rubber compound



Code	Type	Width	Length
AB0010	Bath Mat	400mm	700mm
AB0100	Shower Mat	540mm	540mm

Raised Toilet Seat

- Increases toilet seat height
- Ideal for people with hip or knee injuries or those who find it difficult sitting down or getting up from the toilet
- Easily installed and removed using front location lip and two rear adjustable brackets
- Smooth edges, concealed fittings and high gloss finish make it easier to clean
- Front cut-out for easier personal cleaning
- Moulded plastic is stain and odour resistant



Code	Height	Type	Max User Weight
AF0320	50mm	With Lid	190kg
AF0330	100mm	With Lid	190kg
AF0340	150mm	With Lid	190kg

Ability Commode - AE1890

- Can be used as a bedside commode and overtoilet aid
- Removable toilet seat, lid, commode bowl and back support
- Adjustable seat height
- Plastic moulded armpads for additional support
- Non-slip rubber tips
- Powder coated steel frame
- Plastic commode bowl, toilet seat and lid



Width Between Arms	Seat Height	Arm Height	Base Width	Base Depth	Max User Weight
470mm	410-510mm	170mm	590mm	560mm	125kg

Daily Living

Bathing Sponge - CA0420

- Assists showering or bathing
- Absorbent sponge
- Plastic handle for easy grip
- Long flexible wire stem can be angled
- Eliminates the need for bending and makes reaching easier
- Length 630mm



Stocking Aid - CC0470


- A lightweight flexible plastic gutter to assist in pulling on socks or stockings
- Designed for those who have difficulty in bending or reaching
- Side notches hold the socks or stockings securely
- Sock is placed over gutter then pulled into position over the foot
- Strap ends then pull the sock into place up the leg



Gutter Length	Weight
240mm	45g

Shoe Horn - CC0480

- A long handled shoe horn to assist in putting on shoes
- Designed for those who have difficulty in bending or reaching
- Features a hanging loop attached to the handle
- Supplied with retail packaging



Length	Weight
400mm	45g

Reachers

- Assists in reaching and retrieving a wide variety of household items
- Ideal for those who have difficulty bending or reaching
- Ergonomic full trigger handle requires minimal finger extension
- Ideal for those with weak grip
- Twin gripping jaws are rubber lined for an excellent controlled grip
- Features magnet on head for the retrieval of pins and other small metal objects
- Swivel head allows objects to be picked up from a variety of angles
- Supplied with a plastic clip to attach the reacher to walking frame or stick

Code	Length	Weight
CD0760	600mm	250g
CD0770	750mm	250g
CD0780	1000mm	250g



Handy Soft Turn Pad - LB1970

- Comfortable low profile padded rotating seat
- Reduces the risk of strain or back injury for the carer and user
- Assists the user in transfers from a seated position
- Ideal for transfers in a variety of situations such as car, wheelchair and household areas
- Swivels 360°
- Lightweight and portable
- Cotton upholstered seating surface
- Nylon swivel mechanism provides smooth and safe turning action
- Nonslip underside



Diameter	Height	Weight	Max User Weight
450mm	40mm	350g	125kg

Handy Bar - LB2120

- Provides a safe and supportive handhold for transfers in and out of the car
- Fits securely to cars via the U-shaped striker plate on vehicle door frame
- Contoured handle provides secure grip
- Convenient size and weight
- In-built seat belt cutter and glass breaker for use in an emergency



Handle Length	Handle Width	Weight
210mm	40mm	320g

Pedal Exerciser - HB0040

- Assists greater muscle tone through increased leg movement
- Tension adjustable
- Footstraps prevent feet slipping off pedals



Base Width	Base Depth	Height	Max Pedal Height
400mm	550mm	275mm	320mm

Male Urinal - CF0111

- Autoclavable for sanitising
- Integrated handle
- 100ml graduation markings
- Moulded plastic



Capacity	Weight
1 Litre	115g

Male Urinal with Lid - CF0101

- Integrated handle
- Clip-on lid prevents spillage
- Lid stays attached to urinal while in use
- 100ml graduation markings
- Moulded plastic



Capacity	Weight
1 Litre	165g

Male Urinal Spil-Pruf - CF0141

- Autoclavable for sanitising
- Integrated handle
- Patented non-spill spout prevents spillage in any position
- Spout unscrews easily for disposal and cleaning
- 500ml non-spill capacity
- 100ml graduation markings
- Moulded plastic



Capacity	Weight
1 Litre	240g

Slipper Bed Pan

- Wedge shape assists in positioning pan for use
- Integrated handle
- Front panel prevents spillage when moving or transporting pan
- Autoclavable for sanitising
- Moulded plastic



Code	Size	Width	Length	Height
CF0041	Regular	240mm	320mm	85mm
CF0031	Large	300mm	430mm	90mm

Bed Pan - CF0021

- Frontal opening
- Wide seating area
- Autoclavable for sanitising
- Moulded plastic
- Fits into mobile commode pan carrier



Width	Length	Height
295mm	360mm	115mm

Pressure Care

Apollo Back Support - MB0011

- Medium density foam back support cushion
- Provides comfort and support for the back and spine while sitting
- Maintains correct posture by supporting sideways movement
- Ideal for persons who sit for long periods
- Removable zippered cotton cover



Width	Height
360mm	320mm

All Purpose Cushion - MC0271

- Soft open cell foam
- Convoluted surface for increased air circulation
- Corduroy fabric cover is zippered and removable for cleaning



Width	Depth	Height	Max User Weight
460mm	460mm	60mm	100kg

Ring Cushion, Smooth - MC9630

- Soft open cell foam cushion
- Smooth curved surface
- Ideal for the short term user requiring comfort and support
- Supplied with cotton overslip



Outside Diameter	Inside Diameter	Height	Max User Weight
450mm	130mm	70mm	100kg

Ring Cushion, Convoluted - MC9990

- Soft open cell foam cushion
- Convoluted surface for increased air circulation
- Ideal for the short term user requiring comfort and support
- Supplied with cotton cover



Outside Diameter	Inside Diameter	Height	Max User Weight
430mm	150mm	100mm	100kg


Bed Wedge - BC0680


- Angled foam cushion provides comfort and support to the upper body
- Ideal for those who suffer from respiratory and reflux problems
- Removable navy blue cotton cover
- Wedge can be turned to provide two different recline angles (15° or 75°)
- Foam insert can be removed to alter the angle by approx 5°



Width	Length	Height
680mm	710mm	240 - 600mm







02 9714 1111

General Enquiries - info@carequip.com.au

Orders - orders@carequip.com.au

www.carequip.com.au/smik

5 Millennium Court Silverwater NSW 2128

Proudly owned and distributed by Care Quip

Smik